

# Human CKMT1A/Creatine kinase U-type, mitochondrial ELISA Kit



Catalog No: #EK20009

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CKMT1A/Creatine kinase U-type, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Creatine kinase U-type, mitochondrial
Target Name	Creatine kinase U-type, mitochondrial
Other Names	Acidic-type mitochondrial creatine kinase, Mia-CK, Ubiquitous mitochondrial creatine kinase, U-MtCK, CKMT
Accession No.	P12532
Uniprot	P12532
GeneID	1159;548596;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.401ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CKMT1A

Note: This product is for in vitro research use only