

Human ATF2/Cyclic AMP-dependent transcription factor ATF-2 ELISA Kit



Catalog No: #EK20010

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ATF2/Cyclic AMP-dependent transcription factor ATF-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cyclic AMP-dependent transcription factor ATF-2
Target Name	Cyclic AMP-dependent transcription factor ATF-2
Other Names	cAMP-dependent transcription factor ATF-2, Activating transcription factor 2, Cyclic AMP-responsive element-binding protein 2, CREB-2, HB16, Histone acetyltransferase ATF2, cAMP response element-binding protein CRE-BP1, CREB2, CREBP1
Accession No.	P15336
Uniprot	P15336
GeneID	1386;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.78ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ATF2

Note: This product is for in vitro research use only