

Mouse Atf2/Cyclic AMP-dependent transcription factor ATF-2 ELISA Kit



Catalog No: #EK20011

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Atf2/Cyclic AMP-dependent transcription factor ATF-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclic AMP-dependent transcription factor ATF-2
Target Name	Cyclic AMP-dependent transcription factor ATF-2
Other Names	cAMP-dependent transcription factor ATF-2, Activating transcription factor 2, MXBP protein, cAMP response element-binding protein CRE-BP1
Accession No.	P16951
Uniprot	P16951
GeneID	11909;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Atf2

Note: This product is for in vitro research use only