

Mouse Apoa2/Apolipoprotein A-II ELISA Kit

Catalog No: #EK20012



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Apoa2/Apolipoprotein A-II ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apolipoprotein A-II
Target Name	Apolipoprotein A-II
Other Names	Apo-AII,Apolipoprotein A2
Accession No.	P09813
Uniprot	P09813
GeneID	11807;

Application Details

Assay Type: SandwichDetect Range:1.56-100ng/mL Sensitivity:0.87ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Apoa2

Note: This product is for in vitro research use only