

Human C1QBP/Complement component 1 Q subcomponent-binding protein, mitochondrial ELISA Kit

Catalog No: #EK20013

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

|                    |                                                                                                                                                               |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Human C1QBP/Complement component 1 Q subcomponent-binding protein, mitochondrial ELISA Kit                                                                    |
| Brief Description  | ELISA Kits                                                                                                                                                    |
| Applications       | ELISA                                                                                                                                                         |
| Species Reactivity | Hu                                                                                                                                                            |
| Specificity        | Natural and recombinant Human Complement component 1 Q subcomponent-binding protein, mitochondrial                                                            |
| Target Name        | Complement component 1 Q subcomponent-binding protein, mitochondrial                                                                                          |
| Other Names        | ASF/SF2-associated protein p32,Glycoprotein gC1qBP,C1qBP,Hyaluronan-binding protein 1,Mitochondrial matrix protein p32,gC1q-R protein,p33,GC1QBP,HABP1,SF2P32 |
| Accession No.      | Q07021                                                                                                                                                        |
| Uniprot            | Q07021                                                                                                                                                        |
| GeneID             | 708;                                                                                                                                                          |

Application Details

Assay Type: SandwichDetect Range:78-5000pg/ml Sensitivity:39.9pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:C1QBP

Note: This product is for in vitro research use only