

# Mouse C1qbp/Complement component 1 Q subcomponent-binding protein, mitochondrial ELISA Kit

Catalog No: #EK20014

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse C1qbp/Complement component 1 Q subcomponent-binding protein, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement component 1 Q subcomponent-binding protein, mitochondrial
Target Name	Complement component 1 Q subcomponent-binding protein, mitochondrial
Other Names	GC1q-R protein,Glycoprotein gC1qBP,C1qBP,Gc1qbp
Accession No.	O35658
Uniprot	O35658

## Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:25pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:C1qbp

Note: This product is for in vitro research use only