

Rat C1qbp/Complement component 1 Q subcomponent-binding protein, mitochondrial ELISA Kit

Catalog No: #EK20015

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat C1qbp/Complement component 1 Q subcomponent-binding protein, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Complement component 1 Q subcomponent-binding protein, mitochondrial
Target Name	Complement component 1 Q subcomponent-binding protein, mitochondrial
Other Names	GC1q-R protein,Glycoprotein gC1qBP,C1qBP,Gc1qbp
Accession No.	O35796
Uniprot	O35796
GeneID	29681;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/ml Sensitivity:0.087ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:C1qbp

Note: This product is for in vitro research use only