

Mouse Bche/Cholinesterase ELISA Kit

Catalog No: #EK20016

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Bche/Cholinesterase ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Cholinesterase |
| Target Name | Cholinesterase |
| Other Names | Acylcholine acylhydrolase, Butyrylcholine esterase, Choline esterase II, Pseudocholinesterase |
| Accession No. | Q03311 |
| Uniprot | Q03311 |
| GeneID | 12038; |

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 40.1pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bche

Note: This product is for in vitro research use only