

Pig BCHE/Cholinesterase ELISA Kit

Catalog No: #EK20017



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Pig BCHE/Cholinesterase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Cholinesterase
Target Name	Cholinesterase
Other Names	Acylcholine acylhydrolase,Butyrylcholine esterase,Choline esterase II,Pseudocholinesterase
Accession No.	P32752
Uniprot	P32752

Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.404ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:BCHE

Note: This product is for in vitro research use only