

Rat Agtr1b/Type-1B angiotensin II receptor ELISA Kit

Catalog No: #EK20018

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Agtr1b/Type-1B angiotensin II receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Type-1B angiotensin II receptor
Target Name	Type-1B angiotensin II receptor
Other Names	AT3,Angiotensin II type-1B receptor,AT1B,Agtr1,At1b
Accession No.	P29089
Uniprot	P29089
GeneID	81638;

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.89ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Agtr1b

Note: This product is for in vitro research use only