

Mouse Abo/Histo-blood group ABO system transferase ELISA Kit

Catalog No: #EK20019

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Abo/Histo-blood group ABO system transferase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Histo-blood group ABO system transferase
Target Name	Histo-blood group ABO system transferase
Other Names	Cis-AB transferase,Fucosylglycoprotein 3-alpha-galactosyltransferase,Fucosylglycoprotein alpha-N-acetylgalactosaminyltransferase,Glycoprotein-fucosylgalactoside alpha-N-acetylgalactosaminyltransferase,Glycoprotein-fucosylgalactoside alpha-galactosyltransferase,Histo-blood group A transferase,A transferase,Histo-blood group B transferase,B transferase,NAGAT
Accession No.	P38649
Uniprot	P38649
GeneID	80908;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.161ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Abo

Note: This product is for in vitro research use only