

Mouse Aasdh/Acyl-CoA synthetase family member 4 ELISA Kit



Catalog No: #EK20021

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Aasdh/Acyl-CoA synthetase family member 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Acyl-CoA synthetase family member 4
Target Name	Acyl-CoA synthetase family member 4
Other Names	2-aminoadipic 6-semialdehyde dehydrogenase,Protein LYS2 homolog,Putative aminoadipate-semialdehyde dehydrogenase,Acsf4,U26
Accession No.	Q80WC9
Uniprot	Q80WC9
GeneID	231326;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:40.9pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aasdh

Note: This product is for in vitro research use only