

Mouse Angptl1/Angiopoietin-related protein 1 ELISA Kit

Catalog No: #EK20022

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Angptl1/Angiopoietin-related protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiopoietin-related protein 1
Target Name	Angiopoietin-related protein 1
Other Names	Angiopoietin-like protein 1
Accession No.	Q640P2
Uniprot	Q640P2
GeneID	72713;

Application Details

Assay Type: Sandwich

Detect Range:0.312-20 ng/mL Sensitivity:less than 0.161 ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample

Volume:100ul Assay Time:3-5h Detection wavelength:450 nm

gene_full_name:Angptl1

Note: This product is for in vitro research use only