

Rat Bin1/Myc box-dependent-interacting protein 1 ELISA Kit

Catalog No: #EK20025

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bin1/Myc box-dependent-interacting protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Myc box-dependent-interacting protein 1
Target Name	Myc box-dependent-interacting protein 1
Other Names	Amphiphysin II,Amphiphysin-like protein,Bridging integrator 1,Amph2,Amphl
Accession No.	O08839
Uniprot	O08839
GeneID	117028;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bin1

Note: This product is for in vitro research use only