

## Rat Chia/Acidic mammalian chitinase ELISA Kit

Catalog No: #EK20027

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Chia/Acidic mammalian chitinase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Acidic mammalian chitinase
Target Name	Acidic mammalian chitinase
Other Names	AMCase
Accession No.	Q6RY07
Uniprot	Q6RY07
GeneID	113901;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.075ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Chia

Note: This product is for in vitro research use only