

Mouse Avpr1a/Vasopressin V1a receptor ELISA Kit

Catalog No: #EK20029

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Avpr1a/Vasopressin V1a receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Vasopressin V1a receptor
Target Name	Vasopressin V1a receptor
Other Names	V1aR, AVPR V1a, Antidiuretic hormone receptor 1a, Vascular/hepatic-type arginine vasopressin receptor
Accession No.	Q62463
Uniprot	Q62463
GeneID	54140;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.079ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Avpr1a

Note: This product is for in vitro research use only