

Human ANPEP/Aminopeptidase N ELISA Kit

Catalog No: #EK20033

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Human ANPEP/Aminopeptidase N ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Hu |
| Specificity | Natural and recombinant Human Aminopeptidase N |
| Target Name | Aminopeptidase N |
| Other Names | AP-N,Alanyl aminopeptidase,Aminopeptidase M,AP-M,Microsomal aminopeptidase,Myeloid plasma membrane glycoprotein CD13,gp150,CD13,APN,CD13,PEPN |
| Accession No. | P15144 |
| Uniprot | P15144 |
| GeneID | 290; |

Application Details

Assay Type: Sandwich Detect Range: 7.8-500ng/mL Sensitivity: 4.1ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ANPEP

Note: This product is for in vitro research use only