

## Pig ANPEP/Aminopeptidase N ELISA Kit

Catalog No: #EK20034

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig ANPEP/Aminopeptidase N ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Aminopeptidase N
Target Name	Aminopeptidase N
Other Names	AP-N,Alanyl aminopeptidase,Aminopeptidase M,AP-M,Microsomal aminopeptidase,gp130,CD13
Accession No.	P15145
Uniprot	P15145
GeneID	397520;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.15ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ANPEP

Note: This product is for in vitro research use only