

Mouse C8b/Complement component C8 beta chain ELISA Kit



Catalog No: #EK20035

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse C8b/Complement component C8 beta chain ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Complement component C8 beta chain |
| Target Name | Complement component C8 beta chain |
| Other Names | Complement component 8 subunit beta |
| Accession No. | Q8BH35 |
| Uniprot | Q8BH35 |
| GeneID | 110382; |

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.048ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C8b

Note: This product is for in vitro research use only