

Rat C8b/Complement component C8 beta chain ELISA Kit



Catalog No: #EK20036

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat C8b/Complement component C8 beta chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Complement component C8 beta chain
Target Name	Complement component C8 beta chain
Other Names	Complement component 8 subunit beta
Accession No.	P55314
Uniprot	P55314
GeneID	313421;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.074ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:C8b

Note: This product is for in vitro research use only