

Human B3GAT1/Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1 ELISA Kit



Catalog No: #EK20037

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human B3GAT1/Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1
Target Name	Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1
Other Names	Beta-1,3-glucuronyltransferase 1, Glucuronosyltransferase P, GlcAT-P, UDP-GlcUA:glycoprotein beta-1,3-glucuronyltransferase, GlcUAT-P, GLCATP
Accession No.	Q9P2W7
Uniprot	Q9P2W7
GeneID	27087;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.082ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: B3GAT1

Note: This product is for in vitro research use only