

Rat Actg2/Actin, gamma-enteric smooth muscle ELISA Kit



Catalog No: #EK20038

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Actg2/Actin, gamma-enteric smooth muscle ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Actin, gamma-enteric smooth muscle
Target Name	Actin, gamma-enteric smooth muscle
Other Names	Alpha-actin-3,Gamma-2-actin,Smooth muscle gamma-actin,Acta3,Actsg
Accession No.	P63269
Uniprot	P63269
GeneID	25365;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.123ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Actg2

Note: This product is for in vitro research use only