

Human ADCY9/Adenylate cyclase type 9 ELISA Kit

Catalog No: #EK20039

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ADCY9/Adenylate cyclase type 9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Adenylate cyclase type 9
Target Name	Adenylate cyclase type 9
Other Names	ATP pyrophosphate-lyase 9, Adenylate cyclase type IX, ACIX, Adenylyl cyclase 9, AC9, KIAA0520
Accession No.	O60503
Uniprot	O60503
GeneID	115;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.172ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ADCY9

Note: This product is for in vitro research use only