

Human CSF2RB/Cytokine receptor common subunit beta ELISA Kit



Catalog No: #EK20041

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CSF2RB/Cytokine receptor common subunit beta ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cytokine receptor common subunit beta
Target Name	Cytokine receptor common subunit beta
Other Names	CDw131,GM-CSF/IL-3/IL-5 receptor common beta subunit,CD131,IL3RB,IL5RB
Accession No.	P32927
Uniprot	P32927
GeneID	1439;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.094 ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CSF2RB

Note: This product is for in vitro research use only