

## Mouse Cd63/CD63 antigen ELISA Kit

Catalog No: #EK20042

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cd63/CD63 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CD63 antigen
Target Name	CD63 antigen
Other Names	CD63
Accession No.	P41731
Uniprot	P41731
GeneID	12512;

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cd63

Note: This product is for in vitro research use only