

Rat Bcam/Basal cell adhesion molecule ELISA Kit

Catalog No: #EK20046

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bcam/Basal cell adhesion molecule ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Basal cell adhesion molecule
Target Name	Basal cell adhesion molecule
Other Names	B-CAM cell surface glycoprotein,Lutheran antigen,CD239,Lu
Accession No.	Q9ESS6
Uniprot	Q9ESS6
GeneID	78958;

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:8.2pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bcam

Note: This product is for in vitro research use only