

Mouse Akt2/RAC-beta serine/threonine-protein kinase ELISA Kit



Catalog No: #EK20047

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Akt2/RAC-beta serine/threonine-protein kinase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse RAC-beta serine/threonine-protein kinase
Target Name	RAC-beta serine/threonine-protein kinase
Other Names	Protein kinase Akt-2,Protein kinase B beta,PKB beta,RAC-PK-beta
Accession No.	Q60823
Uniprot	Q60823
GeneID	11652;

Application Details

Assay Type: SandwichDetect Range:15.6-1000pg/mL Sensitivity:7.92pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Akt2

Note: This product is for in vitro research use only