

Mouse Cdkn1a/Cyclin-dependent kinase inhibitor 1 ELISA Kit



Catalog No: #EK20048

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cdkn1a/Cyclin-dependent kinase inhibitor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent kinase inhibitor 1
Target Name	Cyclin-dependent kinase inhibitor 1
Other Names	CDK-interacting protein 1, Melanoma differentiation-associated protein, p21, Cip1, Waf1
Accession No.	P39689
Uniprot	P39689
GeneID	12575;

Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/mL Sensitivity: 37 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: Cdkn1a

Note: This product is for in vitro research use only