

Human CD83/CD83 antigen ELISA Kit

Catalog No: #EK20051

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CD83/CD83 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human CD83 antigen
Target Name	CD83 antigen
Other Names	hCD83,B-cell activation protein,Cell surface protein HB15,CD83
Accession No.	Q01151
Uniprot	Q01151
GeneID	9308;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:40.1pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CD83

Note: This product is for in vitro research use only