

Human B4GALT1/Beta-1,4-galactosyltransferase 1 ELISA Kit



Catalog No: #EK20052

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human B4GALT1/Beta-1,4-galactosyltransferase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Beta-1,4-galactosyltransferase 1
Target Name	Beta-1,4-galactosyltransferase 1
Other Names	Beta-1,4-GalTase 1,UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 1,UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 1,GGTB2
Accession No.	P15291
Uniprot	P15291
GeneID	2683;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:28pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:B4GALT1

Note: This product is for in vitro research use only