

# Mouse B4galt1/Beta-1,4-galactosyltransferase 1 ELISA Kit



Catalog No: #EK20053

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse B4galt1/Beta-1,4-galactosyltransferase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta-1,4-galactosyltransferase 1
Target Name	Beta-1,4-galactosyltransferase 1
Other Names	Beta-1,4-GalTase 1,UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 1,UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 1,Ggtb,Ggtb2
Accession No.	P15535
Uniprot	P15535
GeneID	14595;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.049ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:B4galt1

Note: This product is for in vitro research use only