

# Rat Ccdc82/Coiled-coil domain-containing protein 82 ELISA Kit



Catalog No: #EK20055

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Ccdc82/Coiled-coil domain-containing protein 82 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Coiled-coil domain-containing protein 82
Target Name	Coiled-coil domain-containing protein 82
Accession No.	Q66H73
Uniprot	Q66H73
GeneID	300359;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.182ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Ccdc82

Note: This product is for in vitro research use only