

# Rat Crry/Complement regulatory protein Crry ELISA Kit

Catalog No: #EK20056

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Crry/Complement regulatory protein Crry ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Complement regulatory protein Crry
Target Name	Complement regulatory protein Crry
Other Names	Antigen 5I2
Accession No.	Q63135
Uniprot	Q63135
GeneID	54243;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.183ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Crry

Note: This product is for in vitro research use only