

# Human ASGR2/Asialoglycoprotein receptor 2 ELISA Kit



Catalog No: #EK20057

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ASGR2/Asialoglycoprotein receptor 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Asialoglycoprotein receptor 2
Target Name	Asialoglycoprotein receptor 2
Other Names	ASGP-R 2,C-type lectin domain family 4 member H2,Hepatic lectin H2,HL-2,CLEC4H2
Accession No.	P07307
Uniprot	P07307
GeneID	433;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ASGR2

Note: This product is for in vitro research use only