

Human BMP5/Bone morphogenetic protein 5 ELISA Kit

Catalog No: #EK20058

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BMP5/Bone morphogenetic protein 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Bone morphogenetic protein 5
Target Name	Bone morphogenetic protein 5
Other Names	BMP-5
Accession No.	P22003
Uniprot	P22003
GeneID	653;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BMP5

Note: This product is for in vitro research use only