

Mouse Bmp5/Bone morphogenetic protein 5 ELISA Kit

Catalog No: #EK20059



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Bmp5/Bone morphogenetic protein 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bone morphogenetic protein 5
Target Name	Bone morphogenetic protein 5
Other Names	BMP-5,Bmp-5
Accession No.	P49003
Uniprot	P49003
GeneID	12160;

Application Details

Assay Type: SandwichDetect Range:1.56-100ng/mL Sensitivity:0.78ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Bmp5

Note: This product is for in vitro research use only