

## Human ACVR2A/Activin receptor type-2A ELISA Kit

Catalog No: #EK20060

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ACVR2A/Activin receptor type-2A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Activin receptor type-2A
Target Name	Activin receptor type-2A
Other Names	Activin receptor type IIA,ACTR-IIA,ACVR2
Accession No.	P27037
Uniprot	P27037
GeneID	92;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.087ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ACVR2A

Note: This product is for in vitro research use only