

Rat Acvr2a/Activin receptor type-2A ELISA Kit

Catalog No: #EK20061

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Acvr2a/Activin receptor type-2A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Activin receptor type-2A
Target Name	Activin receptor type-2A
Other Names	Activin receptor type IIA,ACTR-IIA,Actrii,Acvr2
Accession No.	P38444
Uniprot	P38444
GeneID	29263;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Acvr2a

Note: This product is for in vitro research use only