

Human COL4A4/Collagen alpha-4(IV) chain ELISA Kit

Catalog No: #EK20062

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human COL4A4/Collagen alpha-4(IV) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Collagen alpha-4(IV) chain
Target Name	Collagen alpha-4(IV) chain
Accession No.	P53420
Uniprot	P53420
GeneID	1286;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: COL4A4

Note: This product is for in vitro research use only