

Mouse Adc/Antizyme inhibitor 2 ELISA Kit

Catalog No: #EK20063

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adc/Antizyme inhibitor 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Antizyme inhibitor 2
Target Name	Antizyme inhibitor 2
Other Names	AZI2, Arginine decarboxylase-like protein, ADC, ARGDC
Accession No.	Q8BVM4
Uniprot	Q8BVM4
GeneID	242669;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.068ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adc

Note: This product is for in vitro research use only