

## Rat Axin2/Axin-2 ELISA Kit

Catalog No: #EK20065

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Axin2/Axin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Axin-2
Target Name	Axin-2
Other Names	Axin-like protein,Axil,Axis inhibition protein 2,Conductin
Accession No.	O70240
Uniprot	O70240
GeneID	29134;

## Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.312ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Axin2

Note: This product is for in vitro research use only