

Human CH25H/Cholesterol 25-hydroxylase ELISA Kit

Catalog No: #EK20066

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CH25H/Cholesterol 25-hydroxylase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cholesterol 25-hydroxylase
Target Name	Cholesterol 25-hydroxylase
Other Names	Cholesterol 25-monooxygenase,h25OH
Accession No.	O95992
Uniprot	O95992
GeneID	9023;

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.314ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CH25H

Note: This product is for in vitro research use only