

Human ACRV1/Acrosomal protein SP-10 ELISA Kit

Catalog No: #EK20067

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ACRV1/Acrosomal protein SP-10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Acrosomal protein SP-10
Target Name	Acrosomal protein SP-10
Other Names	Acrosomal vesicle protein 1
Accession No.	P26436
Uniprot	P26436
GeneID	56;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ACRV1

Note: This product is for in vitro research use only