

## Rat Afm/Afamin ELISA Kit

Catalog No: #EK20069

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Afm/Afamin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Afamin
Target Name	Afamin
Other Names	Alpha-albumin,Alpha-Alb
Accession No.	P36953
Uniprot	P36953
GeneID	282708;

## Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.38ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Afm

Note: This product is for in vitro research use only