

Mouse Agxt/Serine--pyruvate aminotransferase, mitochondrial ELISA Kit



Catalog No: #EK20070

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Agxt/Serine--pyruvate aminotransferase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Serine--pyruvate aminotransferase, mitochondrial
Target Name	Serine--pyruvate aminotransferase, mitochondrial
Other Names	SPT,Alanine--glyoxylate aminotransferase,AGT,Agxt1
Accession No.	O35423
Uniprot	O35423
GeneID	11611;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Agxt

Note: This product is for in vitro research use only