

Rat Agxt/Serine--pyruvate aminotransferase, mitochondrial ELISA Kit



Catalog No: #EK20071

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Agxt/Serine--pyruvate aminotransferase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Serine--pyruvate aminotransferase, mitochondrial
Target Name	Serine--pyruvate aminotransferase, mitochondrial
Other Names	SPT,Alanine--glyoxylate aminotransferase,AGT,Agt1
Accession No.	P09139
Uniprot	P09139
GeneID	24792;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Agxt

Note: This product is for in vitro research use only