

Human ANTXR1/Anthrax toxin receptor 1 ELISA Kit

Catalog No: #EK20072

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ANTXR1/Anthrax toxin receptor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Anthrax toxin receptor 1
Target Name	Anthrax toxin receptor 1
Other Names	Tumor endothelial marker 8,ATR,TEM8
Accession No.	Q9H6X2
Uniprot	Q9H6X2
GeneID	84168;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ANTXR1

Note: This product is for in vitro research use only