

CENPS Antibody

Catalog No: #44715

Package Size: #44715-1 50ul #44715-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CENPS Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	CENPS antibody detects endogenous levels of total CENPS
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CENPS
Target Name	CENPS
Other Names	APITD1 antibody Apoptosis-inducing TAF9-like domain-containing protein 1 antibody apoptosis-inducing, TAF9-like domain 1 antibody CENP-S antibody CENPS antibody CENPS_HUMAN antibody Centromere protein S antibody FAAP16 antibody FANCM-interacting histone fold protein 1 antibody Fanconi anemia-associated
Accession No.	Swiss-Prot#: Q8N2Z9 NCBI Gene ID:100526739/378708
Uniprot	Q8N2Z9
GeneID	100526739;378708;
Calculated MW	16 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only