

# Human CHRNA7/Neuronal acetylcholine receptor subunit alpha-7 ELISA Kit



Catalog No: #EK20073

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CHRNA7/Neuronal acetylcholine receptor subunit alpha-7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Neuronal acetylcholine receptor subunit alpha-7
Target Name	Neuronal acetylcholine receptor subunit alpha-7
Other Names	NACHRA7
Accession No.	AAH37571/ncbi
Uniprot	AAH375

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.083ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CHRNA7

Note: This product is for in vitro research use only