

Mouse Chrna7/Neuronal acetylcholine receptor subunit alpha-7 ELISA Kit



Catalog No: #EK20074

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Chrna7/Neuronal acetylcholine receptor subunit alpha-7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Neuronal acetylcholine receptor subunit alpha-7
Target Name	Neuronal acetylcholine receptor subunit alpha-7
Other Names	Acra7
Accession No.	P49582
Uniprot	P49582
GeneID	11441;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.39ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Chrna7

Note: This product is for in vitro research use only